



FACE[®] MONOSACCHARIDE COMPOSITION GELS

Product Code: GK60120

This product is intended as a replacement for the sub-components in the FACE Carbohydrate Analysis Kits. Full details for use of this product are included in the protocol supplied with the complete kit.

Supplied Reagents

GK60100 Monosaccharide Composition Gels (10 cm x 10 cm polyacrylamide electrophoresis gels), 10 each

WS-0039 Profiling Gel Reagent Pack

WS-0037 Monosaccharide Running Buffer (formerly Glyko catalog # 70110), 5 each

WS-0040 Mono 2x Loading Buffer (G1) (formerly Glyko catalog # 50169)

WS-0006 Tracking Dye (E2) (formerly Glyko catalog # 50066)

Storage: Shipped ambient for next day delivery. Store only at room temperature upon arrival.

Store reconstituted Monosaccharide Running Buffer at 4°C; may be used for 6 months after reconstitution. Microbial growth in reconstituted buffer is unlikely when stored as indicated. However, since the buffer is not sterile, eventual microbial growth is a possibility. If such growth is observed, discard the buffer and reconstitute a fresh pack.

LICENSE TO USE

Purchase of FACE Carbohydrate Analysis Kits or FACE replacement components from ProZyme or its authorized distributors grants a Use Sublicense under Glyko's FACE patents. By accepting delivery of FACE Kits or replacement components [Material(s)] and by subsequently using them in FACE Carbohydrate Analysis, Recipient agrees to be bound by the following terms and restrictions:

1. A Use Sublicense is granted Recipient for in-house use of Material(s) only.
2. The Material(s) will not be made available by Recipient to any outside parties in any form, separately or in combinations, for any monetary or other consideration or at no charge, except that the Materials may be made available to outside parties who agree to be bound by all the terms and restrictions of this Agreement for purposes of evaluation only.
3. Recipient will not make commercial use of the Material(s) unless it first secures a license agreement from ProZyme, Inc. for such commercial use.
4. Recipient is solely responsible for qualification of the products for the Recipient's specific use.
5. The Material(s) will not be used *in vivo* in humans.

These suggestions and data are based on information we believe to be reliable. They are offered in good faith, but without guarantee, as conditions and methods of use of our products are beyond our control. We recommend that the prospective user determine the suitability of our materials and suggestions before adopting them on a commercial scale.

Suggestions for use of our products or the inclusion of descriptive material from patents and the citation of specific patents in this publication should not be understood as recommending the use of our products in violation of any patent or as permission to license to use any patents of ProZyme, Inc.



1933 Davis Street, Suite 207
San Leandro, CA 94577-1258

TOLL FREE (800) 457-9444
PHONE (510) 638-6900
FAX (510) 638-6919

E-MAIL info@prozyme.com
WEB www.prozyme.com